Neue Therapieoption: PPM-Fisteladapter™

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A new therapeutic option – PPM-fistula adapter™

Introduction of negative pressure wound therapy has improved open abdomen treatment in many ways. Unfortunately handling of intestinal fistulae remains a serious problem. Suction, applied via sponge to the wound surface, often leads to enlargement of the fistula or eventration of the mucosa. Additionally, viscous stool might obstruct the sponge. This results in separation of sponge and wound surface. The suction becomes insufficient and may lead to contamination or infection of the wound.

Several drain-out systems using negative pressure wound therapy have been published. Unfortunately none of them is reliable. Individual, time-consuming modifications are necessary in most cases.

This led to extensive research in the Department of General, Visceral and Vascular Surgery of the University Hospital Magdeburg. The result was the development of the fistula adapter by Dr. Tautenhahn, Prof. Lippert and Primed Halberstadt Medizintechnik GmbH. By now, there are 4 types of the fistula adapter available, for managing enteroatmospheric fistulae.

The original adapter with a height of 3 cm has an inner diameter of 1.5 cm. The lower brim transfers the negative pressure to the wound preventing dislocation. The upper brim serves for fixation of an ostomy bag. The height of the adapter is adjusted to the usual thickness of the negative pressure wound therapy sponge.

In an initial trial 15 patients were treated using the novel fistula adapter. The device was efficient and feasible. Duration for wound dressing change could be reduced considerably. Patients reported improved comfort. Particularly, stability of the wound dressing and possibility of patient mobilization were beneficial. In this trial problems occurred with two of the patients. This led to development of different modifications of the fistula adapter.

For treatment of deep fistulae a 6 cm high adapter with an inner diameter of 1.5cm is now available. For fixation two super-imposed sponges are used. For management of large sole fistulae or multiple small fistulae being close to each other an adapter of 3 cm height and an inner diameter of 4.5 cm can be used. Due to its flexibility it might be fixed in an oval shape.