

**PPM-Wound drainage products**  
**Invasive medical products class IIa**  
**For application with negative pressure from 0 to 500 mmHg**

Description	PZN	Art.-No.	Size	pcs. per unit
<b>PPM drainage port</b> with connecting hose, sterile	08926458	12836	L = 1500 mm	10 x 1
<b>PPM wound foam</b> , anthracite, coarsely porous, PU, sterile	07349409	12801	250 x 150 x 10 mm	10 x 1
	07348083	12802	250 x 150 x 30 mm	10 x 1
	07348108	12803	250 x 300 x 10 mm	10 x 1
	07348114	12804	250 x 300 x 30 mm	10 x 1
	07349148	12809	100 x 150 x 10 mm	10 x 1
	07349378	12810	100 x 150 x 30 mm	10 x 1
<b>PPM abdominal wound foam</b> , anthracite, coarsely porous, PU, sterile	09154644	12813	380 x 250 x 30 mm	1 x 2
<b>PPM wound film</b> , transparent, semipermeable PU film, sterile	10327191	12821	215 x 250 mm	10 x 1
<b>PPM Y-connector</b>	09154650	12815		10 x 1
<b>PPM connector AP</b> , universal adapter for negative pressure generators, large-lock male, sterile	07592713	12827		10 x 1
<b>PPM-Fisteladapter™ Set plus 3/15</b> sterile, composed of: - PPM fistula adapter, 2 pcs. - PPM adapter wound foam, anthracite, coarsely porous - PPM wound film, 3 pcs. - PPM silicone net, 2 pcs.	08929712	12832	L=30 mm, D=15 mm 250 x 300 x 30 mm  215 x 250 mm 200 x 300 mm	1 x 1
<b>PPM-Fisteladapter™ Set plus 3/30</b> sterile, composed of: - PPM fistula adapter, 2 pcs. - PPM adapter wound foam, anthracite, coarsely porous - PPM wound film, 3 pcs. - PPM silicone net, 2 pcs.	01333181	12842	L=30 mm, D=30 mm 250 x 300 x 30 mm  215 x 250 mm 200 x 300 mm	1 x 1
<b>PPM-Fisteladapter™ Set plus 3/45</b> sterile, composed of: - PPM fistula adapter, 2 pcs. - PPM adapter wound foam, anthracite, coarsely porous - PPM wound film, 3 pcs. - PPM silicone net, 2 pcs.	08929729	12833	L=30 mm, D=45 mm 250 x 300 x 30 mm  215 x 250 mm 200 x 300 mm	1 x 1
<b>PPM-Fisteladapter™ Set plus 6/15</b> sterile, composed of: - PPM fistula adapter, 2 pcs. - PPM adapter wound foam, anthracite, coarsely porous, 2 pcs. - PPM wound film, 3 pcs. - PPM silicone net, 2 pcs.	08930069	12834	L=60 mm, D=15 mm 250 x 300 x 30 mm  215 x 250 mm 200 x 300 mm	1 x 1

Description	PZN	Art.-No.	Size	pcs. per unit
<b>PPM-Fisteladapter™ Set standard 3/15</b> <b>sterile</b> , composed of: - PPM-fistula adapter, 2 pcs. - PPM-adapter-wound foam, anthracite, coarsely porous - PPM wound film , 3 pcs.	7592742	12825	L=30 mm, D=15 mm 250 x 300 x 30 mm  215 x 250 mm	1 x 1
<b>PPM-Fisteladapter™ Set standard 3/30</b> <b>sterile</b> , composed of: - PPM fistula adapter, 2 pcs. - PPM adapter wound foam, anthracite, coarsely porous - PPM wound film, 3 pcs.	1333152	12841	L=30 mm, D=30 mm 250 x 300 x 30 mm  215 x 250 mm	1 x 1
<b>PPM-Fisteladapter™ Set standard 3/45</b> <b>sterile</b> , composed of: - PPM fistula adapter, 2 pcs. - PPM adapter wound foam, anthracite, coarsely porous - PPM wound film, 3 pcs.	7592759	12826	L=30 mm, D=45 mm 250 x 300 x 30 mm  215 x 250 mm	1 x 1
<b>PPM-Fisteladapter™ Set standard 6/15</b> <b>sterile</b> , composed of: - PPM fistula adapter, 2 pcs. - PPM adapter wound foam, anthracite, coarsely porous, 2 pcs. - PPM wound film, 3 pcs.	7751956	12829	L=60 mm, D=15 mm 250 x 300 x 30 mm  215 x 250 mm	1 x 1
<b>PPM silicone net, sterile</b> transparent silicone coated wound dressing	8930075	12835	200 x 300 mm	1 x 1

Last update 01/2015